


# iQ-DASHBOARD

PRESENTS

YOUR iQ-WEB PACS STATISTICS AT A GLANCE\*



SELECT FILTERS AND TIME RANGES.

A BAR CHART WILL THEN PROVIDE YOUR STATISTICS.

## FILTER OPTIONS:

- Referring physician
- Modality
- Modality and body part
- Body part and modality
- Body part
- Site
- Site and modality
- AE title
- Duplicate patient IDs
- Studies without SR
- Studies without PDF

## TIME RANGES:

- Today
- Last night
- Yesterday
- Two nights ago
- Last 7 days
- Last 30 days
- This month
- Previous month
- Last 12 months
- This year
- Previous year
- FROM-TO

# ADVANTAGES

- ▶ Easy performance reviews:
  - ▶ For PACS administrators to timely detect any technical issues
  - ▶ For financial or radiology directors to evaluate productivity and service quality
- ▶ Data export to CSV files for the creation of advanced statistics using Microsoft Excel
- ▶ User management to give access to selected users
- ▶ Individual filter preferences can be saved as default by each user

## LIVE MONITOR

...shows all incoming studies of the last 30 minutes.

Each study includes a link to the respective study in the PACS, so the PACS administrator can make any required organizational changes, such as forwarding, deleting or moving the study.

SITE	STUDY DATE	DATE OF BIRTH	PATIENT ID	ACCESSION NUMBER	PATIENT NAME
Wayne and partners radiology	2016-10-31 12:13:57	1958-05-24	<a href="#">P.200645</a>	ACC_P.200645S.2	Lastname_200645 Middelname_200645 Firstname_200645
St. Jones clinic	2016-10-31 00:07:02	1998-01-19	<a href="#">P.200707</a>	ACC_P.200707S.0	Lastname_200707 Middelname_200707 Firstname_200707
St. Jones clinic	2016-10-29 03:06:57	1998-01-19	<a href="#">P.200710</a>	ACC_P.200710S.0	Lastname_200710 Middelname_200710 Firstname_200710
Sacred Heart	2016-10-29 00:07:23	1955-07-16	<a href="#">P.200712</a>	ACC_P.200712S.0	Lastname_200712 Middelname_200712 Firstname_200712
Wayne and partners radiology	2016-10-30 17:10:11	1960-01-24	<a href="#">P.200747</a>	ACC_P.200747S.1	Lastname_200747 Middelname_200747 Firstname_200747
Sacred Heart	2016-10-30 02:55:48	1960-01-24	<a href="#">P.200747</a>	ACC_P.200747S.0	Lastname_200747 Middelname_200747 Firstname_200747
Sacred Heart	2016-10-28 20:04:01	1987-12-03	<a href="#">P.200746</a>	ACC_P.200746S.0	Lastname_200746 Middelname_200746 Firstname_200746
Memorial State University	2016-11-02 10:26:28	1978-06-04	<a href="#">P.200745</a>	ACC_P.200745S.1	Lastname_200745 Middelname_200745 Firstname_200745